



This tutorial is provided solely as a procedure workflow example. Solutions, quantities, containers and equipment should be adjusted or replaced to suit the needs of your application(s) and each specimen processed.

IMPORTANT

Before proceeding with the following Steps, remove reagents from the refrigerator to ensure that the media, used to reduce non-specific binding and retain cell morphology, and paramagnetic antibodies are at room temperature.

Step 1

(Optional)

Prior to target cell enrichment, prepare a slide smear from the bone marrow specimen for morphological review. The slide smear will serve as a reference that sufficient viable cells are present.

Use an appropriate cellular stain to identify the presence of target cells on your smear.

Step 2

Using a 15 mL conical tube, create a homogeneous mixture by combining RPMI 1640 media with the bone marrow specimen.

- Depending on plasma cell percentage, allocate 100-200 μ L of bone marrow for each EpiSep[®] HS (<3%=200 μ L, \geq 3%=100 μ L).
- As specified in the chart below, use RPMI 1640 media to bring the specimen/media mixture volume to a ratio of 1 mL for each EpiSep HS.

Specimen/media volume to EpiSep HS conversion chart:

1 mL = 1 EpiSep HS	6 mL = 6 EpiSep HS
2 mL = 2 EpiSep HS	7 mL = 7 EpiSep HS
3 mL = 3 EpiSep HS	8 mL = 8 EpiSep HS
4 mL = 4 EpiSep HS	9 mL = 9 EpiSep HS
5 mL = 5 EpiSep HS	10 mL = 10 EpiSep HS

Step 3

- Vortex to resuspend the paramagnetic antibodies.
- As specified in the chart below, add 10 μ L of the appropriate paramagnetic antibody to the specimen/media mixture in the 15 mL tube for each EpiSep HS.

Paramagnetic antibody volume to EpiSep HS conversion chart:

10 μ L = 1 EpiSep HS	60 μ L = 6 EpiSep HS
20 μ L = 2 EpiSep HS	70 μ L = 7 EpiSep HS
30 μ L = 3 EpiSep HS	80 μ L = 8 EpiSep HS
40 μ L = 4 EpiSep HS	90 μ L = 9 EpiSep HS
50 μ L = 5 EpiSep HS	100 μ L = 10 EpiSep HS

Step 4

Cap and gently invert the 15 mL conical tube several times (3-4). Attach the 15 mL CellCycler Tube Adapter to the 15 mL conical tube then insert it into the CellCycler. Set voltage to 7.5V DC (14.0 - 16.0 RPM) for 15-30 minutes. The CellCycler is used to keep cells and reagents in uniform suspension during the incubation of cells with paramagnetic antibodies.

Step 5

- At the end of the incubation period, loosen the cap of the 15 mL conical tube and place it firmly into a MTD-15 (#A4102-1) magnetic tube dock for 5 minutes. **NOTE:** Do not incubate more than 5 minutes in the MTD-15 as it may induce aggregation.
- At the end of the 5-minute incubation, remove supernatant from the side opposite the MTD-15 magnet. Be careful not to scrape cells off the side of the conical tube.
- After the supernatant is removed with an aspirator or transfer pipette, place the supernatant into a waste container labeled biohazard. **NOTE:** Supernatant contains non-target cells that may be useful for further analysis.
- As specified below, resuspend target cells with 5-10 mL of hypotonic solution (0.075M KCl). Cap the 15 mL conical tube, invert to mix cells and place in a 37°C environment for 10 minutes.
Specimen volumes < 500 μ L = 5 mL of 0.075M KCl.
Specimen volumes \geq 500 μ L = 10 mL of 0.075M KCl.

Step 6

- At the end of the 10-minute incubation period, remove the 15 mL conical tube from the incubator. Slowly add 1 mL of Carnoy's fixative. Invert gently to mix.
- Loosen the cap of the 15 mL conical tube and place it firmly into a MTD-15 (#A4102-1) magnetic tube dock for 5 minutes. **NOTE:** Do not incubate more than 5 minutes in the MTD-15 as it may induce aggregation.
- At the end of the 5-minute incubation, remove supernatant from the side opposite the MTD-15 magnet. Be careful not to scrape cells off the side of the conical tube.
- Resuspend the cells with Carnoy's fixative using the same volume conversion in Step #2. **NOTE:** For storage of the tube for reflex testing, please refer to www.wavesense.com/faq

Step 7

Label each EpiSep HS with a pencil and insert into the MSD5 (#A1104-1) five position magnetic slide dock. Make sure it clicks to indicate full insertion.

Step 8

- Add 1 mL of suspension from the 15 mL conical tube to the inner wall of each EpiSep HS. Allow 1 minute for the suspension to drain or until all fluid has drained from the well. **NOTE:** Do not let the slides sit and dry with the EpiSep HS cap still attached.
- Remove EpiSep HS cap by gently pulling up while slide is still wet and in the MSD5. Place slides into a humidity-controlled environment for 5 minutes. (45-50% humidity)

Step 9

- Slides are now ready for your laboratory's standard FISH protocol or special stain process. Inspect monolayered cell target area by viewing slides under a phase contrast microscope.